

# Identify<sup>®</sup> CUP TESTS

\* NOTE: Panel intended for forensic use only.

ID-CP5 ID-CP6 ID-H6-1 ID-H7-1 ID-H9-1 ID-CP10 ID-CP12-BUP ID-CP12-BUP-ADULT ID-US12-BUP MADE IN USA ID-H12-EFKK-1 ID-H12-1 ID-H12-2 ID-H13-1 ID-US14-ADULT MADE IN USA ID-H14-ADULT-1 ID-H14-ADULT-EFKT-2 ID-H16-ADULT-EFKT

SUBSTANCE TESTED	CUP TESTS																
	ID-CP5	ID-CP6	ID-H6-1	ID-H7-1	ID-H9-1	ID-CP10	ID-CP12-BUP	ID-CP12-BUP-ADULT	ID-US12-BUP	ID-H12-EFKK-1	ID-H12-1	ID-H12-2	ID-H13-1	ID-US14-ADULT	ID-H14-ADULT-1	ID-H14-ADULT-EFKT-2	ID-H16-ADULT-EFKT
Amphetamines (AMP)		1000	500	500	500	1000	1000	1000	500	500	500	500	500	500	500	500	500
Barbiturates (BAR)					300	300	300	300	300		300	300	300	300	300		300
Benzodiazepines (BZO)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Buprenorphine (BUP)				10	10		10	10	10	10	10	10	10	10	10	10	10
Cocaine (COC)	300	300	150	150	150	300	300	300	150	150	150	150	150	150	150	150	150
Ecstasy (MDMA)							500	500	500		500	500	500	500	500	500	500
Marijuana (THC)	50	50	50	50		50	50	50	50	50	50	50	50	50	50	50	50
Methadone (MTD)					300	300	300	300	300		300	300	300	300	300	300	
Methadone Metabolite (EDDP)														300	300		300
Methamphetamines (MET)					500	1000	1000	1000	500	500	500	500	500	500	500	500	500
Opiates (MOP) or (OPI300)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Oxycodone (OXY)	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100
Phencyclidine (PCP)						25	25	25	25		25		25	25	25		25
Tricyclic Antidepressants (TCA)												1000	1000	1000	1000		
* Ethyl Glucuronide (ETG)										500						500	500
* Fentanyl (FEN)										25						20	25
* Kratom (KRA)										100							
* Synthetic Marijuana (K2)										25						25	25
* Tramadol (TRA)																100	100

## ADULTERANTS

Creatinine (C)											✓					✓	✓	✓	✓
Glutaraldehyde (G)																✓	✓		✓
Acidic or Alkaline pH (P)											✓					✓	✓	✓	✓
Nitrite (N)																✓	✓		✓
Oxidants (O)																✓	✓	✓	✓
Specific Gravity (S)											✓					✓	✓		✓

12 PANEL



Identify<sup>®</sup> Diagnostics

14 PANEL



Identify<sup>®</sup> Diagnostics USA

14 PANEL



Identify<sup>®</sup> Health


# Identify<sup>®</sup> DIP TESTS

\* NOTE: Panel intended for forensic use only.


ID-H1-COT-DIP  
ID-H1-ETG-DIP  
ID-H1-FEN-DIP  
ID-H1-GAB-DIP  
ID-H1-KRA-DIP  
ID-CP5-DIP  
ID-CP10-DIP  
ID-H10-DIP  
ID-CP12-DIP  
ID-H12-EFT-DIP

SUBSTANCE TESTED	DIP TESTS							
Amphetamines (AMP)				1000	1000	500	1000	500
Barbiturates (BAR)					300	300	300	
Benzodiazepines (BZO)					300	300	300	300
Buprenorphine (BUP)						10	10	10
Cocaine (COC)				300	300	150	300	150
Ecstasy (MDMA)					500		500	
Marijuana (THC)				50	50	50	50	50
Methadone (MTD)					300	300	300	300
Methamphetamines (MET)				1000	1000	500	1000	500
Opiates (MOP) or (OPI300)				300	300	300	300	300
Oxycodone (OXY)					100	100	100	100
Phencyclidine (PCP)							25	
* Cotinine / Nicotine (COT)	200							
* Ethyl Glucuronide (ETG)		500						500
* Fentanyl (FEN)			25					50
* Gabapentin (GAB)				2000				
* Kratom (KRA)				100				
* Tramadol (TRA)								100

5 PANEL  
single sided



12 PANEL  
double sided



1 PANEL  
(FEN)  
fentanyl



1 PANEL  
(ETG)  
alcohol



1 PANEL  
(GAB)  
gabapentin



12 PANEL  
double sided

