



DRUGS TESTED

ALL CUP TESTS with cut-off levels

	ID-CP5	ID-CP6	ID-US6	ID-US6-BUP	ID-H16-1	ID-H19-1	ID-CP10	ID-H11-1	ID-CP12-BUP	ID-CP12-BUP-ADULT	ID-US12-BUP	ID-H12-1	ID-CP13	ID-US14-ADULT	ID-H14-ADULT-EFKT	ID-H16-ADULT-EFKT
DRUGS TESTED	CUP TESTS															
Amphetamines (AMP)		1000	500		500	500	1000	500	1000	1000	500	500	1000	500	1000	500
Barbiturates (BAR)						300	300	300	300	300	300	300	300	300	300	300
Benzodiazepines (BZO)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Buprenorphine (BUP)				10		10		10	10	10	10	10	10	10	10	10
Cocaine (COC)	300	300	150	150	150	150	300	150	300	300	150	150	300	150	300	150
Ecstasy (MDMA)								500	500	500	500	500	500	500	500	500
Marijuana (THC)	50	50	50		50		50		50	50	50	50	50	50	50	50
Methadone (MTD)				300		300	300	300	300	300	300	300	300	300	300	300
Methadone Metabolite (EDDP)														300		300
Methamphetamines (MET)						500	1000	500	1000	1000	500	500	1000	500	1000	500
Opiates (MOP) or (OPI)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	2000	300
Oxycodone (OXY)	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100
Phencyclidine (PCP)							25	25	25	25	25	25	25	25		25
Tricyclic Antidepressants (TCA)												1000	1000			
** Cotinine / Nicotine (COT)																
** Ethyl Glucuronide (ETG)															500	500
** Fentanyl (FEN)															25	25
** Synthetic Marijuana (K2)															25	25
** Tramadol (TRA)															100	100
*ADULTERANTS																
Creatinine (C)										✓				✓	✓	✓
Glutaraldehyde (G)														✓		✓
Acidic or Alkaline pH (P)										✓				✓	✓	✓
Nitrite (N)														✓		✓
Oxidants (O)														✓	✓	✓
Specific Gravity (S)										✓				✓		✓

** NOTE: Premium drug panel not CLIA Waived intended for forensic use only. * NOTE: Drug tests with adulterants measure levels of substances in human urine used to tamper or falsify a urine test result.

A green box  indicates the drug or adulterant in the far left column is tested for. The number inside the box indicates the cut-off level in ng/ml. A cut-off level is the concentration a drug must meet or exceed in the urine sample to trigger a positive result. Learn more about cut-off levels & detection times later in this image gallery.

Tests above with the American Flag  are designed & manufactured in the United States. All Identify Diagnostics & Identify Diagnostics USA tests are made with USA reagents. Reagents are the substances on each panel strip which produce the drug test result.

DRUGS TESTED

ALL DIP TESTS with cut-off levels

	ID-COT-DIP	ID-THC-DIP	ID-H1-DIP-THC	ID-H1-DIP-COC	ID-H1-DIP-ETG	ID-H1-DIP-FEN	ID-CP5-DIP	ID-CP6-DIP	ID-CP7-DIP-ADULT	ID-H8-DIP-1	ID-CP10-DIP	ID-US10-DIP-1	ID-CP12-DIP
DRUGS TESTED	DIP TESTS												
Amphetamines (AMP)						1000			500	1000	500	1000	
Barbiturates (BAR)										300	300	300	
Benzodiazepines (BZO)							300	300	300	300	300	300	
Buprenorphine (BUP)									10	10		10	10
Cocaine (COC)				150		300	300	300	150	300	150	300	
Ecstasy (MDMA)											500		500
Marijuana (THC)		50	50			50	50	50		50	50	50	
Methadone (MTD)										300	300	300	300
Methadone Metabolite (EDDP)													
Methamphetamines (MET)						1000	1000	1000	500	1000	500	1000	
Opiates (MOP) or (OPI)						300	300	300	300	300	300	300	
Oxycodone (OXY)							100	100	100	100	100	100	
Phencyclidine (PCP)													25
Tricyclic Antidepressants (TCA)													
** Cotinine / Nicotine (COT)	200												
** Ethyl Glucuronide (ETG)				500									
** Fentanyl (FEN)					25								
** Synthetic Marijuana (K2)													
** Tramadol (TRA)													
*ADULTERANTS													
Creatinine (C)									✓				
Glutaraldehyde (G)													
Acidic or Alkaline pH (P)										✓			
Nitrite (N)													
Oxidants (O)										✓			
Specific Gravity (S)													

HEP-C



HIV



LC-MS SUPPLY



CORFLEX BRACES



CareFull
catch
BY CAREFULL PRODUCTS



Target Markets

In Office Physician Labs- CLIA
Employment Testing
Parental / Family Testing
Government
Addiction

Distributed By

Medical Distribution Group, Inc.
2361 Whitfield Park Avenue
Sarasota, FL 34243
P: 727-744-2967
F: 727-499-7355